dopey dictionaries

2024.11.25

The data from prodigy is saved to a .xy file which is then loaded and sorted by the dopey.load() method for visualization and initial analysis.

Meta data

Meta data keys that are present in (almost) all dopey dicts:

file_name	string	The name of the loaded file.
spectrum_id	integer	An id set by Prodigy.
experiment	dict	Contains information about the measurement, such as
		lens mode, etc. See below.
type	string	Type of measurement as identified by dopey.
labels	dict	Labels for axes and intensity data. See below.

The experiment dict

This is an example of the experiment dict for a Fermi map measurement, i.e. a deflector map.

```
{'Version': '4.99.1-r110443',
 'Energy_Axis': 'Kinetic energy (eV)',
 'Count Rate': 'Intensity (counts/s)',
 'Spectrum ID': 13,
 'Analysis Method': 'UPS',
 'Analyzer': 'PhoibosCCD',
 'Lens Mode': 'WideAngleMode',
 'Scan Mode': 'SnapshotFAT',
 'Curves Per Scan': 200,
 'Values Per Curve': 1,
 'Dwell Time': 0.25,
 'Excitation Energy': 24.0,
 'Ek': 7.29821,
 'Ep': 7.0,
 'Ordinate Range': [-15.0, 15.0],
 'parameters': ['ShiftX [a.u.]'],
 'Column_labels': ['energy', 'counts/s']}
```

The labels dict

This is an example of the labels dict for the same measurement as above:

```
{'x': 'Kinetic energy (eV)',
  'y': 'Detector y-range',
  'z': 'X-Deflector (deg.)',
  'intensity': 'Intensity (counts/s)'}
```

Data

Data keys includes axes and intensity. E.g., an ARPES cut has two axes (x and y, where x is energy and y is analyzer angles) and a 2d intensity array, while a Fermi map has three axes (x, y, and z, where z is the deflector angle) and a 3d intensity array.

X	array	Axis (1d)	
У	array	Axis (1d) if applicable	
Z	array	Axis (1d) if applicable	
intensity	string	Array (1d, 2d, 3d)	

Data types

dopey will try to identify what kind of measurement and data that is contained in the file when loaded and then sort the data accordingly. At the moment the following types of data can be identified (continuously updated):

<i>Type</i> xps	Axes X	energy	<i>Intensity</i> [y, x]
	У	y-direction	
arpes	x y	energy analyzer angle	[y, x]
fermi_map	x y z	energy analyzer angle deflector angle	[z, y, x]
target_scattering_spectrum	X	energy	[x]
spin_edc (*)	x polarity	energy target magn. dir.	[p, scan, x]
spin_mdc (*)	x y polarity	energy deflector angle target magn. dir.	[p, scan, y, x]
spin_map (*)	x y z polarity	energy deflector angle 1 deflector angle 2 target magn. dir.	[p, scan, z, y, x]
spin_arpes	x y	energy deflector angle	[y, x]

For the three spin measurements spin_edc, spin_mdc, and spin_map there is an additional intensity key intensity_mean. This contains an intensity array with one dimension lower than intensity.

Dicts

xps

file_name string spectrum_id integer experiment dict type string Χ array (1d) energy array (1d) detector y-range intensity array (2d) labels dict

arpes

file_name string spectrum_id integer experiment dict type string array (1d) Χ energy array (1d) detector angle У intensity array (2d) labels dict

fermi_map

file_name string spectrum_id integer dict experiment type string Χ array (1d) energy detector angle array (1d) У deflector angle Z array (1d) intensity array (3d) labels dict

spin_edc

file_name string
spectrum_id integer
experiment dict
type string
x array (1d) energy

polarity array (1d)
intensity array (3d)
intensity_mean array (2d)
labels dict

spin_mdc

file_name string spectrum_id integer experiment dict string

x array (1d)

energy

energy

detector angle deflector angle

detector angle

y array (1d) polarity array (1d)

intensity array (1d) array (4d) array (3d) labels dict

spin_map

file_name string spectrum_id integer experiment dict type string

x array (1d)

y array (1d) z array (1d)

polarity array (1d)
intensity array (5d)
intensity_mean array (4d)
labels dict

spin_arpes

file_name string spectrum_id integer experiment dict string

x array (1d) energy

y array (1d) deflector angle

intensity array (4d) labels dict